



Breast Ultrasound Optimization eLearning

The Breast Ultrasound Optimization eLearning course will provide an overview of breast ultrasound image optimization. Topics covered include a review of breast anatomy, image quality and optimization technical features available to improve the ultrasound image quality. The course will also cover some of the more common pitfall encountered in breast ultrasound.

The Breast Ultrasound Optimization eLearning module will:

- Describe the normal anatomy of the breast and the sonographic appearance
- Review the sonographic terminology of breast ultrasound
- Recognize the various optimization technical features available in breast ultrasound imaging
- Identify the more common pitfalls of breast ultrasound
- Associate benefits of great image quality with diagnostic confidence and workflow

This course is intended for customers who wish to acquire knowledge and skill in performing breast ultrasound.

This ASRT-accredited module is a self-paced tutorial consisting of enlightening content and relevant images.



For more information

Contact a Philips ultrasound clinical services coordinator at 800-522-7022 and visit our education catalog at

www.learningconnection.philips.com/ultrasound