



**PHILIPS**

Clinical Education

# Philips Ultrasound DICOM<sup>®</sup> Encoding

Our goal at Philips Healthcare is to provide the clinical education you need to make the most of your equipment investment.

This tutorial discusses the Philips Ultrasound DICOM Encoding application. The module reviews the history of DICOM and its relevance in imaging.

The course will describe the tool used to create custom DICOM elements which can be transferred to PACS.

Upon completion of this module, you will be able to:

- Describe DICOM
- Explain what DICOM enables
- Describe a Data Dictionary
- Describe a DICOM Conformance Statement
- Summarize the concept of DICOM encoding
- Explain the user-defined measurements / calculations creation and transfer processes

#### **For more information**

Contact Philips Ultrasound Clinical Education at 800.522.7022 and visit our education catalog at [www.learningconnection.philips.com/ultrasound](http://www.learningconnection.philips.com/ultrasound)

Please visit [www.learningconnection.philips.com/ultrasound](http://www.learningconnection.philips.com/ultrasound)

**PHILIPS**