



Philips EPIQ Ultrasound Basic System Training Curriculums

Our goal at Philips Healthcare is to provide the clinical education you need to make the most of your equipment investment.

The aim of this comprehensive series is to provide the physician, sonographer, and echo cardiographers with a comprehensive library of courses at an exceptional value.

Purchase the educational series directly from the Philips Learning Center by clicking Register Today next to each curriculum.

We are pleased to offer three NEW e-learning curriculums focusing on the Philips EPIQ Ultrasound System.

In today's healthcare climate, it is critical to use your ultrasound system to its fullest potential. Philips courses are designed to educate customers on the leading edge technologies specific to each system and application.

These ASRT-accredited modules are self-paced tutorials consisting of enlightening content and colorful images. There are no prerequisites for these courses. You may register for an individual course or register for an entire curriculum.

For questions about classes or purchasing call: 1-800-522-7022 or visit www.learningconnection.philips.com/ultrasound

[Philips EPIQ Ultrasound Basic System Training Curriculum \(3.5 CEUs\)](#)

- System Overview
 - System Setups and Annotations
 - Patient Data Management and Connectivity
 - Advanced Support Service
- REGISTER TODAY**

[Philips EPIQ Ultrasound Basic System Training Curriculum for General Imaging and Vascular \(4 CEUs\)](#)

- Analysis Setups
 - Spectral Doppler
 - 2D Optimization
 - Color Doppler
 - Using Analysis
- REGISTER TODAY**

[Philips EPIQ Ultrasound Basic System Training Curriculum for Cardiovascular \(4.25 CEUs\)](#)

- Analysis Setups
 - Contrast Imaging
 - Stress Echo Protocols
 - Tissue Doppler Imaging
 - 2D Optimization
 - Color and Spectral Doppler Imaging
 - Using Analysis
- REGISTER TODAY**