



PHILIPS

Clinical Education

Ultrasound Lab Accreditation Quality Control

Our goal at Philips Healthcare is to provide the clinical education you need to make the most of your equipment investment.

This tutorial discusses ultrasound lab accreditation requirements with a focus on the quality control component. Learners will review the components of a continuous quality control program and see examples of proper documentation.

The Ultrasound Lab Accreditation Quality Control module will:

- Introduce the three major organizations that offer ultrasound lab accreditation
- Review Quality Control (QC) program guidelines
- Explain proper documentation of quality control tests
- Discuss phantom testing procedures

For more information

Contact Philips Ultrasound Clinical Education at 800.522.7022 and visit our education catalog at www.learningconnection.philips.com/ultrasound

Please visit www.learningconnection.philips.com/ultrasound