



PHILIPS

Clinical Education

Philips X8-2t TEE Transducer eLearning

This eLearning course provides an overview of the unique features and benefits of the Philips x8-2t transducer.

Clinical tips, an image optimization review and helpful comments from experienced and knowledgeable TEE course instructors are also included in this review.

This eLearning module will:

- Review Philips History of Innovation in TEE
- Discuss the features of the X8-2t
- Explain useful tips on using the X8-2t
- Review image optimization in Live 3D
- Explain changes to viewing and cropping modes
- Review tips from our Live 3D TEE course instructors

This ASRT-accredited module is a self-paced tutorial consisting of enlightening content and colorful images.

For more information

Contact Philips Ultrasound Clinical Education at 800.522.7022 and visit our education catalog at www.learningconnection.philips.com/ultrasound

Please visit www.learningconnection.philips.com/ultrasound