



PHILIPS

Clinical Education

TEE Electrical Leakage Testing

We know it is critical that your ultrasound transducers operate efficiently and that your sonographers are properly trained to today's standards. Ultrasound laboratory accreditation requires the structural and electrical integrity of the TEE transducer must be checked between each use, using an ultrasound transducer leakage tester. Our Philips Clinical Education Specialists value the opportunity to train your staff to perform the electrical leakage testing. The benefit for you is better patient care.

The TEE Electrical Leakage Testing course covers the reasons for performing leakage testing and the recommended steps in the process.

Our Philips Clinical Specialist may demonstrate the use of the customer's ULT - 2000 testing device during the course.

Course Objectives

- Discuss the rationale behind transducer leakage testing
- Explain the inspection process for TEE transducers
- Review the steps for performing a leakage test
- Discuss troubleshooting the testing process

For more information

Contact Philips Ultrasound Clinical Education at 800.522.7022 and visit our education catalog at www.learningconnection.philips.com/ultrasound

Please visit www.learningconnection.philips.com/ultrasound

© 2016 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. (Royal Philips) or their respective owners.

PHILIPS